



100%
Made in Italy
since 1955

Labò 2.0 EC PF



1

Filtro
in uscita
HEPA

HEPA
exhaust
filter



2

Filtro a
cartuccia
classe M

Cartridge
filter, class M
certified



3

Pompa
automatica
scuoli-filtro

Automatic
filter
cleaning



Unità di aspirazione mobile adatta all'impiego nel laboratorio odontotecnico per l'aspirazione delle polveri di lavorazione e dei residui abrasivi. Collegabile a sabbiatrici, lucidatrici, fresatrici, micromotori, ed utilizzabile con i banchi da laboratorio.
@ Sistema di filtraggio a 3 stadi (secco in microfibra, filtro a cartuccia classe M, filtro uscita HEPA)
@ Elettronica EMC e presa elettroutensile con sistema di avviamento e spegnimento automatico
@ Pulizia automatica e manuale del filtro

Mobile extraction unit suitable for use in the dental laboratory for the suction of dust, abrasives and residues. Connectable to sand-blasters, polishers, milling machines, micromotors and for use in the laboratory bench.
@ Filtering 3 stages system (micropor bag, Class M cartridge filter, HEPA exhaust filter)
@ EMC control and power-tool socket with automatic start and stop system
@ Automatic and manual filter cleaning

Modello Model	Motore Motor	Depressione Depression	Rumorosità Noise level	Dimensioni Dimension	Peso Weight
EC- PF	900 W	1900 mm/H ₂ O	59 dBA max	24 x 31, h 48 cm	15 kg

Your complete range of dust collectors